

### 3. Before Applying for After Service

#### The remote control fails to turn on?

Q : Is the power switch on? A: Turn on the power switch.

Q: Is the power cord properly connected? A: Connect the power cord.

Q: Is the remote control properly connected? A: Connect the remote control.

Q: Is the fuse out? A: Replace the fuse.

#### The external projector fails to operate?

Q: Are the External Projectors properly connected? A: Connect external the projector.

Q: IAre the External Projectors hot? A: If the temperature display of the remote control indicates 'tH,'"disconnect the power cord and contact your dealer or service center (Do not disassemble the product).

#### The product is malfunctioning?

Q: Is the remote control properly connected? A: Connect the remote control.

Q: Is the temperature display of the remote indicating SE, tH, CE, LE, od? A: Check the remote control error message on page 17 and contact your dealer or service center.

Q: Is the internal projector not responding to the ↑ ↓ buttons? A: The ↑ ↓ buttons only function in semi-auto mode. No response in semi-auto mode requires you to contact a service center.

Q: When turned on, does the Internal Projector continuously move unup and down? A: Disconnect the power cord and contact your dealer or a service center.